

Concentric Rings water baths

The **SESS** Refrigerated Circulating baths with a special control module provide excellent accuracy & uniformity essential for most demanding applications in the temperature range of -10°C to + 120° C. These heating cooling baths are ideal for both immersing samples & circulating controlled temperature liquid to external equipments.

Applications:

Concentric Rings water baths are identical for all applications requiring gentle steaming work with Erlenmeyer flasks, glasses etc.

Advantages:

- Fitted with a water level indicating device
- Optional constant water level device.
- A removable flat drop on lid for concentric rings enables the bath to be used for other work also.
- Double walled, insulated, full stainless steel construction.
- Concentric rings in steps of 25 and 50 mm ID moulded from industrial grade non-corrosive material.
- A precise Digital Temperature Indicator cum Controller having resolution and control accuracy of $\pm 1^\circ\text{C}$ in digital model.



Model	Bath	Control	Concentric Rings	Heater
CR 6	30 x 20 x 10 cm	Digital	6 holes of 75 mm Ø with steps of 25 mm	800 W x 1
CR12	40 x 30 x 10 cm	Digital	12 holes of 75 mm Ø with steps of 25 mm	800 W x 2
CR 24	60 x 40 x 10 cm	Digital	24 holes of 75 mm Ø with steps of 25 mm	1500 W x 2

SRI EASWARI SCIENTIFIC SOLUTION

Old No- 499, New no: 2/20, Plot No - 110A, Srinivasanagar, 2nd cross street, Srinivasapuram, Baraniputhure, Iyyapanthangal,
Near: srinivasa welding shop, Chennai-600 056 INDIA,