

Platform shaker

The **SESS** Platform shaker incorporate design features that offer maintenance-free shaking under a constant and reproducible temperature environment. After preparing your samples from growing cultures that express proteins to staining your gels for visualization, trust our orbital shaker for reliable around-the-clock operations to keep your cells alive and flourishing within your working environment.

Applications:

Shaker are ideal for all applications requiring precise and reproducible results for Incubations, fermentations, homogenization, chemical and biochemical reactions, enzyme and tissue studies

Advantages:

- Digital RPM as standard
- Constant shaking speed even during heater operation.
- Double walled, insulated, full stainless steel construction.
- Shaking Speed: up to 300 RPM.



Technical Data: Specification / Model No	SESS-PS
Lotus Clamps	16 x 250ml (or) 16 x 500ml for conical flask
Tray size	420 x 420 mm
Shaking Tray	Stainless Steel 304 Grade
Outer	Mild Steel Powder coated
Timer	Provided with Digital Timer
RPM Indicator	Provided
Shaking Speed	50 – 300 RPM - ORBITAL
Motor	PMDC Motor
Shaking dia	25mm
Power Supply	230V, 50HZ, Single Phase

SRI EASWARI SCIENTIFIC SOLUTION

Old No- 499, New no: 2/20, Plot No - 110A, Srinivasanagar, 2nd cross street, Srinivasapuram, Baraniputhure, Iyyapanthangal,
Near: srinivasa welding shop, Chennai-600 056 INDIA,