

## Unstirred Water Bath

The **SESS** unstirred bath are offered in digital as well as thermostatic control. These small compact unstirred water baths are constructed in full stainless steel & are provided with special energy efficient heaters.

### *Applications:*

Ideal for sample pre-treatment in QC, for Serological tests in Pathology and all thermo stating tests

### *Advantages:*

- Fully stainless steel, including bottom, makes it GMP standard & total rust proof.
- Easy to set & use - direct digital setting & readout.
- Double walled full SS construction with glass wool insulation.
- Temperature Range: Ambient + 5 to 65°C.
- Digital Temperature Controller cum Indicator - 1°C resolution.
- Energy efficient tubular elmech heater.



Model	Bath	Control	Heater
UWB40	27 x 13 x 12 cm	Digital	300 W x 1
UWB100	30 x 27 x 12 cm	Digital	500 W x 1

### SRI EASWARI SCIENTIFIC SOLUTION

Old No- 499, New no: 2/20, Plot No - 110A, Srinivasanagar, 2<sup>nd</sup> cross street, Srinivasapuram, Baraniputhure, Iyyapanthangal,  
Near: srinivasa welding shop, Chennai-600 056 INDIA,