

## HOT AIR OVEN - UNIVERSAL DIGITAL

- Temperature range : 50°C to 250°C
- Temperature Accuracy :  $\pm 1^\circ\text{C}$
- Temperature Controller : Digital Microprocessor controller with LED Display
- Heating : Side Heated
- Air Circulation : Blower
- Inner chamber : SS 304
- Cabinet : GI powder coated
- Outer Door Gasket : Silicone
- Trays : Perforated SS Tray
- Supply : 220-230 volts AC



Model No.	Internal Dim. W x D x H cm	Capacity Ltrs	No of Trays	Watts
HO-01UD	30 x 30 x 30	27	1	1000
HO-02UD	35 x 35 x 35	42	2	1200
HO-03UD	45 x 45 x 45	91	2	1500
HO-04UD	45 x 45 x 60	121	3	1800
HO-05UD	60 x 60 x 60	216	3	2250
HO-06UD	60 x 60 x 90	324	3	2250

### SRI EASWARI SCIENTIFIC SOLUTION PVT. LTD.

Door No.2/298, ANE Garden Permul Kovil Street, Srinivasapuram  
Paraniputhur Post, Iyyappanthangal, Chennai - 600 122, Tamilnadu, India.



ISO 9001 : 2015

SL.NO	OTHER TERMS AND CONDITIONS:	
1	PAYMENT:	70 % Advance along with purchase Order and balance 30% on After completion of Inspection at our factory against Proforma invoice.
2	OUR BANK DETAILS:	Name: Indian Overseas Bank, Branch: Iyyappanthangal, Chennai. A/C number: 228002000027748, IFSC Code: IOBA0002280, MICR Code: 600020176.
3	GST	
	SALES:	Currently 18%. (Extra as applicable at the time of dispatch or billing).
	SERVICE:	Currently 18%. (Extra as applicable at the time of dispatch or billing). Installation charges extra.
4	WARRANTY:	Against any manufacture defect - 12 months from the date of Installation and Commissioning or 13 months from date of dispatch of material, whichever is earlier?
5	PRICE BASIS:	Ex-Works, Chennai, India.
6	PACKING AND FORWARDING:	Extra in this price.
7	FRIEGHT:	Extra in this price.
8	INSURANCE:	Not in our scope. To be arranged by you, if required based on our dispatch details.
9	DELIVERY:	8-10 Weeks from the date of receipt of Techno - commercially confirmed order with advance payment.
10	VALIDITY:	60 days
11	JURISDICTION:	This quotation along with Terms & Conditions mentioned herein is subject to Chennai Jurisdiction only.
12	EXCLUSION:	Masonry and carpentry works, civil works, electrical and hydraulic connections, disassembling of pre-existent installations, transport from unloading point to installation place, lifting devices, materials such as oxygen, acetylene etc. taxes & duties on the supply.
13	ORDER AMENDMENTS:	Any changes in the specifications which amounts to change in model and Qty shall entail revision of delivery & prices
14	FORCE MAJEURE:	No liability shall be attached to us for non-performance or delayed execution of the order as a result of "force majeure"
15	LOSS OF DAMAGE IN TRANSIST:	Our responsibility ceases on delivering the goods to the carriers and we undertake no liabilities for any damage or loss in transit whatsoever. All claims for shortage or damage must reach us within 48 hours in case of local supply or within 15 days if it is an upcountry supply. Post delays excludes after these periods have elapsed. No claim will be recognized by us.
16	NOTE:	The product undergoes continuous R&D and up gradation and Design specifications are subject to change without notice.

For SRI EASWARI SCIENTIFIC SOLUTION PVT LTD,

B.PRIYA  
Mobile : 9444427728|  
E: [info@sess.co.in](mailto:info@sess.co.in)

**SRI EASWARI SCIENTIFIC SOLUTION PVT. LTD.**

Door No.2/298, ANE Garden Permul Kovil Street, Srinivasapuram  
Paraniputhur Post, Iyyappanthangal, Chennai - 600 122, Tamilnadu, India.



ISO 9001 : 2015